

MaxAbsorb RAPID

Polyglycolic Acid Fast Suture

Properties

- Retention of 50% knot tensile strength after 6 to 8 days
- Complete mass absorption within 42 days

Indications

- Gynaecology/Obstetric procedures
- Urology procedures
- Skin closure

Advantages

- High initial tensile strength
- Predictable and constant degradation rate
- Good knot security
- Excellent handling properties
- Quick mass absorption



Product Overview

Presentation/Format	: Needled suture packed in a Dozen box
Structure	: Multifilament, Braided
Colour	: Undyed
Chemical composition	: Polyglycolic acid
Sizes	: USP 5/0 (metric 1) - USP 1 (metric 4)
Type of absorption	: Hydrolysis
Sterilization	: Gamma radiation