

Maximus Surgical Sutures - Veterinary Surgery



Contents

Wound Closure

MaxiMus Veterinary Surgery Needles

Veterinary Suture Types

Understanding Veterinary Suture Codes

Decoding the codes

Material Prefixes

Suffixes

Veterinary Needle Specifications

Veterinary Needle Shapes

Abbreviation for specification

Needle Point Symbol

MAXIMUS & MAXIMUS CAS

MAXABSORB & MAXABSORB CAS

MAXIMUS & MAXIMUS PLAIN CAS

MAXIMUS CHROMIC & MAXIMUS CHROMIC CAS

MAXIMIDE & MAXIMIDE CAS

MAXPROP & MAXPROP CAS







MaxiMus Veterinary Surgery Needles

In pursuit of offering best quality products, Our sutures are paired with the superior stainless steel needles for an unparalleled surgical experience.

Our sutures come with sharper and stronger needles stainless steel specifically designed to meet the ever growing needs in all the surgical disciplines

With every needle swaged using Multi Stroke Technology (MST) and tested for a needle pull check, Our needles are distinguished with a high penetration coating and a square body design for a firm grip.



Needle Profiles	Point type	Symbol
	Taper point	
	Conventional cutting	
	Reverse cutting	
	Round bodied/Taper point	
	Straight	

Veterinary Suture Types


Suture Type	Size	Usage	Support Duration
Absorbable Sutures			
Polyglactin 910	1, 0, 2-0, 3-0, 4-0	used to suture fascia, muscles, or ligature of blood vessels.	21-28 days
Polyglycolic Acid Sutures	1, 0, 2-0, 3-0, 4-0	used to suture fascia, muscles, or ligature of blood vessels.	21-28 days
Catgut Sutures	1, 0, 2-0, 3-0, 4-0	used to suture mucous membrane.	07-10 days
Chromic Catgut Sutures	1, 0, 2-0, 3-0, 4-0	used to suture fascia, muscles, or ligature of blood vessels	10-20 days

Suture Type	Size	Usage
Non Absorbable Sutures		
Silk	3-0, 4-0, 5-0	used to suture mucous membrane
Nylon	3-0, 4-0, 5-0	most commonly used to close and suture skin after surgery or trauma to the skin.
Polypropylene	4-0,5-0,6-0	used to suture nerve, tendon or blood vessels.

Understanding Veterinary Suture Codes

Our code logic represents the necessary information needed to identify the required suture-needle combination. The suture codes are alphanumeric in nature to indicate the suture material, the suture size, and the reference code. The alphabetical suffixes are a key to understand the additional features about the suture and/or the needle.

Decoding the Codes




















Material Prefixes

ABBREVIATION	MEANING
PGN	Polyglactin 910
PGNMC	Polyglactin 910 Cassette Suture
PGAMC	Polyglycolic Acid Cassette Suture
PCL	Poliglecaprone 25
PDX	Polydioxanone
PLGMC	Plain Catgut Cassette Suture
CGL	Chromic Catgut Loop
CHGMC	Chromic Catgut Cassette Suture
SLK	Silk
SLS	Silk Non Needled
NYL	Polyamide
NYS	Polyamide Non Needled
NYLMC	Nylon (Polyamide) Cassette Suture
PPL	Polypropylene
PPG	Polypropylene with Pledgets
PPLMC	Polypropylene Cassette Suture


Suffixes

ABBREVIATION	MEANING
VB	Black Needle
BX	Black Needle; Extra Length Suture
DN	Double Needle
DNL	Double Needle; Long Length Suture
DS	Doctor Special Code
G	Green Suture
ML	Multi Strands
SL	Short Length Suture
U	Undyed Suture
V	Violet Suture
W	White Suture
X	Extra Pair of Suture
XF	Extra Fine Needle
XL	Extra Length Suture

Veterinary Needle Specifications

Type	Taper	Reverse cutting	Diamond	Cutting	Premium cutting	Spatula	Taper Cutting
Design							
Point							
Shape							
1/4 Circle		3/8 Circle			1/2 Circle		

Veterinary Needle Shapes



Abbreviation for Specification

ABBREVIATION	MEANING
HC	1/2 CIRCLE
QC	1/4 CIRCLE
CU	3/8 CIRCLE
FE	5/8 CIRCLE
RC	REVERSE CUTTING
CT	CUTTING
RB	ROUND BODIED
TC	TAPER CUTTING
ST	STRAIGHT
SM	SPATULA
VIO	VIOLET
UND	UNDYED

Needle Point Symbol



MaxiMus CAS

Polyglactin 910 Suture

Properties

- Retention of 75% knot tensile strength after 14 days
- Retention of 50% knot tensile strength after 21 days
- Complete mass absorption within 56 to 70 days

Indications

- Gastrointestinal procedures
- Gynaecology procedures
- Ophthalmic procedures
- Orthopedic procedures
- Vascular Surgery
- Skin closure (intracutaneous, subcutaneous)
- Ligatures
- Connective Tissues
- Muscles

Advantages

- Smooth knot run-down/Knot slide
- Gentle knot repositioning
- Excellent knot security
- Minimal inflammatory reaction



Product Overview

Presentation/Format	: Needled suture packed in a Dozen box Non-needed suture packed in a cassette
Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
Chemical composition	: Poly (glycolide/l-lactide) 90/10
Coating	: Poly(glycolide/l-lactide) + Calcium Stearate
Sizes	: Dyed - USP 8/0 (metric 0.4) - USP 2 (metric 5) : Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilization	: Ethylene oxide




Description	Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 2-0 X 90-55MM ST CT	PGN209713	2-0		55 mm Straight Conventional Cutting	90	12	VIO
PGLA VIO 2-0 X 90-45MM CU RC	PGN209282	2-0		45 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 90-30MM HC RB	PGN20 2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 140-30MM HC RB TC DN	PGN20 2317DN	2-0		30 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGLA VIO 2-0 X 70-30MM HC RB	PGN20 2341	2-0		30 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 90-40MM HC RB	PGN20 2345	2-0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 90-40MM HC RC	PGN20 2382	2-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA UND 2-0 X 90-26MM CU CT	PGN20 2404U	2-0		26 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 2-0 X 45-26MM HC RB	PGN200775	2-0		26 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 75-26MM HC RB	PGN209121	2-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 70-26MM HC RB	PGN209123	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 75-31MM HC RB	PGN209136	2-0		31 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-36MM HC RB	PGN209140	2-0		36 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-40MM HC RB	PGN209150	2-0		40 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-65MM HC RB	PGN209158	2-0		65 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-26MM HC TC	PGN20 9350	2-0		26 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-40MM HC TC	PGN209359	2-0		40 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-31MM HC TC	PGN209360	2-0		31 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-36MM HC TC	PGN209363	2-0		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-45MM HC TC	PGN209375	2-0		45 mm 1/2 Circle Taper Cut	75	12	VIO





Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 2-0 X 90-36MM HC RB	PGN209440	2-0		36 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 90-36MM HC TC	PGN209463	2-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 2-0 X 90-40MM HC TC	PGN209466	2-0		40 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 2-0 X 100-31MM HC RB	PGN209637	2-0		31 mm 1/2 Circle Round Body	100	12	VIO
PGLA UND 2-0 X 75-40MM HC TC	PGN209904	2-0		40 mm 1/2 Circle Taper Cut	75	12	UND
PGLA UND 2-0 X 75-36MM HC RB	PGN200258	2-0		36 mm 1/2 Circle Round Body	75	12	UND
PGLA UND 2-0 X 70-26MM HC RB	PGN200417	2-0		26 mm 1/2 Circle Round Body	70	12	UND
PGLA VIO 2-0 X 70-36MM HC RC	PGN200466	2-0		36 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 2-0 X 90-30MM CU RC	PGN200589	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA UND 2-0 X 70-26MM HC RC	PGN200869	2-0		26 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 2-0 X 45-17MM HC RB	PGN200306	2-0		17 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 70-21MM HC RB	PGN200824	2-0		21 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-22MM HC RC	PGN20461	2-0		22 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 2-0 X 70-24MM CU RC	PGN200443	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 2-0 X 45-36MM HC RC	PGN200710	2-0		36 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA UND 2-0 X 70-36MM HC RB BP	PGN200259	2-0		36 mm 1/2 Circle Blunt Point	70	12	UND
PGLA VIO 2-0 X 70-36MM FE RB	PGN200375	2-0		36 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-40MM HC RC	PGN200478	2-0		40 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 2-0 X 70-40MM HC RB	PGN200351	2-0		40 mm 1/2 Circle Round Body	70	12	VIO

Description	Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 2-0 X 75-30MM HC RB	PGN200368	2-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 75-30MM HC TC	PGN200517	2-0		30 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 75-35MM HC RB BLUNT	PGN200738	2-0		35 mm 1/2 Circle Blunt Point	75	12	VIO
PGLA VIO 2-0 X 90-35MM CU RC	PGN209390	2-0		35 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 75-25MM HC CT	PGN209225	2-0		25 mm 1/2 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-45MM HC CT	PGN209245	2-0		45 mm 1/2 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-35MM CU CT	PGN209235	2-0		35 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-26MM CU CT	PGN209226	2-0		26 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-24MM CU RB	PGN209224	2-0		24 mm 3/8 Circle Round Body	75	12	VIO
PGLA VIO 2-0 X 70-24MM CU RC	PGN200453	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA VIO 2-0 X 90-48MM HC RB	PGN209290	2-0		48 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 2-0 X 90-36MM HC TC	PGN209464	2-0		36 mm 1/2 Circle Taper Cut	90	12	UND
PGLA VIO 2-0 X 70-31MM HC RC	PGN209231	2-0		31 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA UND 2-0 X 90-36MM HC RB	PGN20 945U	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 2-0 X 90-30MM CU RC	PGN20 2316	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 2-0 X 75-22MM HC TC	PGN20 9828	2-0		22 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 2-0 X 45-30MM HC RB	PGN20 2317SL	2-0		30 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 2-0 X 90-26MM HC RB	PGN20 2356	2-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 35-26MM HC RB	PGN20 2356SL	2-0		26 mm 1/2 Circle Round Body	35	12	VIO
PGLA VIO 2-0 X 70-26MM HC RB VB	PGN20 2122	2-0		26 mm 1/2 Circle Round Body	70	12	VIO






Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 2-0 X 90-30MM HC RB DN	PGN20 2318DN	2-0		30 mm 1/2 Circle Round Body Round Body	90	12	VIO
PGLA VIO 2-0 X 45-19MM SKKI RB	PGN20 2903	2-0		19 mm SKKI Round Body	45	12	VIO
PGLA UND 2-0 X 76-60MM ST CT	PGN20 2390U	2-0		60 mm Straight Conventional Cutting	76	12	UND
PGLA VIO 2-0 X 70-36MM HC RB	PGN20 323	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 70-26MM FE RB	PGN20 602	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 90-60MM ST CT	PGN209718	2-0		60 mm Straight Conventional Cutting	90	12	VIO
PGLA UND 2-0 X 75-26MM CU CT	PGN202404USL	2-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA VIO 2-0 X 75-36MM J TYPE RB	PGN209220J	2-0		36 mm J Shape Round Body	75	12	VIO
PGLA UND 2-0 X 75-36MM CU CT	PGN209532	2-0		36 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 2-0 X 75-26MM CU RC	PGN209572	2-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA UND 2-0 X 45-26MM CU CT	PGN209522	2-0		26 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 2-0 X 75-60MM ST CT	PGN2013360	2-0		60 mm Straight Conventional Cutting	75	12	VIO
PGLA VIO 2-0 X 75-35MM CU RC	PGN209390SL	2-0		35 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 2-0 X 90-26MM HC RB EX	PGN200369	2-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 2-0 X 70-26MM FE RB EX	PGN200602	2-0		26 mm 5/8 Circle Round Body	70	12	VIO
PGLA UND 2-0 X 90-36MM HC RB EX	PGN200945	2-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 2-0 X 70-36MM HC RB EX	PGN200339	2-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 2-0 X 75-31MM J Type RB	PGN20 9136J	2-0		31 mm J Shape Round Body	75	12	VIO
PGLA UND 2-0 X 70-51MM ST CT	PGN209273	2-0		51 mm Straight Conventional Cutting	70	12	UND






















Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 2-0 X 75-26MM HC RB-3D	PGN209121Z	2-0		26 mm 1/2 Circle Round Body	75	36	VIO
PGLA UND 2-0 X 70-24MM CU RC-3D	PGN200443Z	2-0		24 mm 3/8 Circle Reverse Cutting	70	36	UND
PGLA VIO 2-0 X 70-36MM HC RC-3D	PGN200466Z	2-0		36 mm 1/2 Circle Reverse Cutting	70	36	VIO
PGLA VIO 2-0 X 75-36MM HC RB-3D	PGN209140Z	2-0		36 mm 1/2 Circle Round Body	75	36	VIO
PGLA VIO 2-0 X 90-22MM HC RB-3D	PGN200312XLZ	2-0		22 mm 1/2 Circle Round Body	90	36	VIO
PGLA UND 2-0 X 70-40MM HC RC	PGN200267U	2-0		40 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA UND 2-0 X 75-55MM ST CT	PGN209713U	2-0		55 mm Straight Conventional Cutting	75	12	UND
PGLA VIO 2-0 X 75-17MM HC RB	PGN200306XL	2-0		17 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 2-0 X 70-22MM HC RC	PGN206934U	2-0		22 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA UND 2-0 X 90-40MM HC RB	PGN200957U	2-0		40 mm 1/2 Circle Round Body	90	12	UND
PGLA UND 3-0 X 90-25MM FE RB	PGN309730	3-0		25 mm 5/8 Circle Round Body	90	12	UND
PGLA VIO 3-0 X 90-25MM CU CT	PGN302643	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	VIO
PGLA UND 3-0 X 90-25MM CU CT	PGN30 2328U	3-0		25 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 3-0 X 45-22MM CU CT	PGN30 2401	3-0		22 mm 3/8 Circle Conventional Cutting	45	12	VIO
PGLA UND 3-0 X 90-19MM CU CT	PGN30 2402U	3-0		19 mm 3/8 Circle Conventional Cutting	90	12	UND
PGLA VIO 3-0 X 70-20MM HC RB	PGN30 2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 45-16MM CU CT	PGN30 2471	3-0		16 mm 3/8 Circle Conventional Cutting	45	12	VIO
PGLA VIO 3-0 X 90-22MM HC CT	PGN30 2472	3-0		22 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGLA VIO 3-0 X 90-36MM HC TC	PGN30 2516	3-0		36 mm 1/2 Circle Taper Cut	90	12	VIO
PGLA VIO 3-0 X 75-26MM HC RB	PGN309120	3-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 75-31MM HC RB	PGN30 9331N	3-0		31 mm 1/2 Circle Round Body	75	12	VIO























Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA UND 3-0 X 45-19MM CU RC	PGN304970	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 3-0 X 90-17MM HC TC	PGN305180	3-0		17 mm 1/2 Circle Taper Cut	90	12	UND
PGLA UND 3-0 X 45-16MM HC RC	PGN300657	3-0		16 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA VIO 3-0 X 70-17MM HC RB	PGN300305	3-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 45-19MM CU CT	PGN300824	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 3-0 X 70-19MM CU RC	PGN300427	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 45-22MM HC RC	PGN3079013	3-0		22 mm 1/2 Circle Reverse Cutting	45	12	VIO
PGLA VIO 3-0 X 70-22MM HC RC	PGN300460	3-0		22 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 3-0 X 70-22MM HC RB	PGN300311	3-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 90-24MM CU RC	PGN300442	3-0		24 mm 3/8 Circle Reverse Cutting	90	12	UND
PGLA VIO 3-0 X 45-26MM CU RC	PGN3078713	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 3-0 X 70-26MM HC RC	PGN300868	3-0		26 mm 1/2 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 90-26MM HC RB	PGN300527	3-0		26 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 3-0 X 90-36MM HC RB	PGN300344	3-0		36 mm 1/2 Circle Round Body	90	12	UND
PGLA VIO 3-0 X 90-40MM HC RC	PGN30 9389	3-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGLA VIO 3-0 X 90-40MM HC RB	PGN30 9388	3-0		40 mm 1/2 Circle Round Body	90	12	VIO
PGLA VIO 3-0 X 75-30MM HC RB	PGN30864	3-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 75-25MM HC RB	PGN300528	3-0		25 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 70-25MM HC RB	PGN300526	3-0		25 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 75-17MM HC TC	PGN300216	3-0		17 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 3-0 X 90-35MM HC RB	PGN300738	3-0		35 mm 1/2 Circle Round Body	90	12	VIO



Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 3-0 X 90-25MM CU RB	PGN309325	3-0		25 mm 3/8 Circle Round Body	90	12	VIO
PGLA VIO 3-0 X 90-26MM HC CT	PGN309326	3-0		26 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGLA VIO 3-0 X 75-26MM CU RC	PGN3078713XL	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 3-0 X 70-60MM ST CT	PGN300664	3-0		60 mm Straight Conventional Cutting	70	12	VIO
PGLA VIO 3-0 X 70-19MM CU RC	PGN300398	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA UND 3-0 X 75-30MM CU RC	PGN309330U	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA VIO 3-0 X 75-30MM CU RC	PGN309329	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGLA VIO 3-0 X 70-30MM HC RC	PGN309330	3-0		30 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGLA VIO 3-0 X 70-26MM HC RB VB	PGN30 2123	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 45-17MM HC TC	PGN30 2515	3-0		17 mm 1/2 Circle Taper Cut	45	12	VIO
PGLA VIO 3-0 X 70-12MM HC RB	PGN30 2427	3-0		12 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 3-0 X 90-26MM HC RB DN	PGN30 2436DN	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	VIO
PGLA VIO 3-0 X 45-12MM CU RC	PGN30 2495	3-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA VIO 3-0 X 70-36MM HC RB	PGN30 322	3-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 45-24MM CU RC	PGN30 683U	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 3-0 X 70-24MM CU RC	PGN30 936U	3-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 45-17MM HC RB	PGN300774	3-0		17 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 3-0 X 90-18MM HC RB	PGN309318	3-0		18 mm 1/2 Circle Round Body	90	12	VIO
PGLA UND 3-0 X 75-60MM ST CT	PGN309717	3-0		60 mm Straight Conventional Cutting	75	12	UND
PGLA VIO 3-0 X 75-20MM HC RB	PGN309114	3-0		20 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 3-0 X 45-22MM HC CT	PGN309472	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	VIO
PGLA UND 3-0 X 45-26MM CU CT	PGN309521	3-0		26 mm 3/8 Circle Conventional Cutting	45	12	UND










Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA UND 3-0 X 75-19MM CU CT	PGN309571	3-0		19 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 3-0 X 75-26MM CU CT	PGN309581	3-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 3-0 X 75-26MM CU RC	PGN309890	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGLA UND 3-0 X 75-17MM HC TC	PGN309826	3-0		17 mm 1/2 Circle Taper Cut	75	12	UND
PGLA VIO 3-0 X 45-16MM CU RC	PGN300608	3-0		16 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA VIO 3-0 X 90-25MM FE RB	PGN309525	3-0		25 mm 5/8 Circle Round Body	90	12	VIO
PGLA UND 3-0 X 45-16MM CU CT	PGN309507	3-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 3-0 X 75-36MM CU CT	PGN309531	3-0		36 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA VIO 3-0 X 75-16MM CU CT	PGN309444	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 3-0 X 70-36MM HC RB EX	PGN300338	3-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 3-0 X 70-24MM CU RC EX	PGN3093605	3-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 3-0 X 70-26MM CU RC	PGN30 2407	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA UND 3-0 X 45-24MM CU RC EX	PGN300683	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 3-0 X 45-22MM HC CT	PGN309741	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	UND
PGLA VIO 3-0 X 75-26MM CU CT	PGN309581V	3-0		26 mm 3/8 Circle Conventional Cutting	75	12	VIO
PGLA VIO 3-0 X 90-26MM HC RB-3D	PGN300527Z	3-0		26 mm 1/2 Circle Round Body	90	36	VIO
PGLA UND 3-0 X 90-24MM CU RC-3D	PGN300442Z	3-0		24 mm 3/8 Circle Reverse Cutting	90	36	UND
PGLA UND 3-0 X 45-19MM CU RC-3D	PGN304970Z	3-0		19 mm 3/8 Circle Reverse Cutting	45	36	UND
PGLA VIO 3-0 X 75-26MM HC TC	PGN309327	3-0		26 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA VIO 3-0 X 70-22MM HC RB-3D	PGN300311Z	3-0		22 mm 1/2 Circle Round Body	70	36	VIO
PGLA VIO 3-0 X 75-26MM HC RC	PGN309375	3-0		26 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGLA UND 3-0 X 70-60MM ST CT	PGN300663U	3-0		60 mm Straight Conventional Cutting	70	12	UND
PGLA VIO 3-0 X 75-20MM HC TC	PGN309118	3-0		20 mm 1/2 Circle Taper Cut	75	12	VIO
PGLA UND 3-0 X 70-30MM CU RC	PGN300502U	3-0		30 mm 3/8 Circle Reverse Cutting	70	12	UND

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------










MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA UND 3-0 X 45-26MM HC RB	PGN300864U	3-0		26 mm 1/2 Circle Round Body	45	12	UND
PGLA VIO 3-0 X 135-36MM HC RB RC DN	PGN309336DN	3-0		36 mm 1/2 Circle Round Body Reverse Cutting	135	12	VIO
PGLA VIO 4-0 X 90-16MM CU RC	PGN409437	4-0		16 mm 3/8 Circle Reverse Cutting	90	12	VIO
PGLA VIO 4-0 X 90-13MM ST CT	PGN409410	4-0		13 mm Straight Conventional Cutting	90	12	VIO
PGLA VIO 4-0 X 70-20MM HC RB	PGN40 2304	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGLA VIO 4-0 X 70-16MM HC RB	PGN40 2305	4-0		16 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 4-0 X 70-16MM CU CT	PGN40 2443U	4-0		16 mm 3/8 Circle Conventional Cutting	70	12	UND
PGLA VIO 4-0 X 45-12MM CU RC	PGN40 2494	4-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 4-0 X 70-19MM CU RC	PGN400422	4-0		19 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA UND 4-0 X 45-19MM CU RC	PGN400496	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 4-0 X 45-08MM HC RC	PGN400504	4-0		8 mm 1/2 Circle Reverse Cutting	45	12	UND
PGLA VIO 4-0 X 45-12MM HC RB	PGN4074313	4-0		12 mm 1/2 Circle Round Body	45	12	VIO
PGLA UND 4-0 X 45-12MM HC RB	PGN400434	4-0		12 mm 1/2 Circle Round Body	45	12	UND
PGLA UND 4-0 X 45-16MM CU CT	PGN400845	4-0		16 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 4-0 X 45-17MM HC RB	PGN400714	4-0		17 mm 1/2 Circle Round Body	45	12	UND
PGLA UND 4-0 X 70-17MM HC RB	PGN400214	4-0		17 mm 1/2 Circle Round Body	70	12	UND
PGLA VIO 4-0 X 45-22MM HC RB	PGN4077113	4-0		22 mm 1/2 Circle Round Body	45	12	VIO
PGLA VIO 4-0 X 70-22MM HC RB	PGN400310	4-0		22 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 4-0 X 45-24MM CU RC	PGN400682	4-0		24 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA UND 4-0 X 70-24MM CU RC	PGN400935	4-0		24 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA UND 4-0 X 70-26MM CU RC	PGN400441	4-0		26 mm 3/8 Circle Reverse Cutting	70	12	UND
PGLA VIO 4-0 X 70-24MM CU RC	PGN400451	4-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA VIO 4-0 X 45-26MM HC RB	PGN400773	4-0		26 mm 1/2 Circle Round Body	45	12	VIO

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 4-0 X 70-26MM HC RB	PGN400315	4-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGLA UND 4-0 X 45-13MM CU RC	PGN400464	4-0		13 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA VIO 4-0 X 75-16MM HC RB	PGN400715	4-0		16 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 4-0 X 75-16MM CU RB	PGN400846	4-0		16 mm 3/8 Circle Round Body	75	12	VIO
PGLA VIO 4-0 X 60-13MM CU RB	PGN400465	4-0		13 mm 3/8 Circle Round Body	60	12	VIO
PGLA VIO 4-0 X 70-19MM CU RC	PGN400397	4-0		19 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGLA VIO 4-0 X 70-17MM HC TC	PGN409417	4-0		17 mm 1/2 Circle Taper Cut	70	12	VIO
PGLA VIO 4-0 X 75-26MM HC RB	PGN409425	4-0		26 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 4-0 X 90-16MM CU CT HCRB DN	PGN40 2354DN	4-0		16 mm 3/8 Circle 1/2 Circle Conventional Cutting Round Body	90	12	VIO
PGLA VIO 4-0 X 90-45MM ST CT	PGN40 2388	4-0		45 mm Straight Conventional Cutting	90	12	VIO
PGLA VIO 4-0 X 45-16MM FE CT	PGN40 2465	4-0		16 mm 5/8 Circle Conventional Cutting	45	12	VIO
PGLA VIO 4-0 X 70-19MM ST RB	PGN40 2573	4-0		19 mm Straight Round Body	70	12	VIO
PGLA VIO 4-0 X 75-17MM HC RB	PGN40 9106	4-0		17 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 4-0 X 90-16MM CU CT	PGN40 2354	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	VIO
PGLA VIO 4-0 X 45-16MM CU RC	PGN409437SL	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA VIO 4-0 X 75-20MM HC RB	PGN409113	4-0		20 mm 1/2 Circle Round Body	75	12	VIO
PGLA VIO 4-0 X 45-13MM ST CT	PGN409410SL	4-0		13 mm Straight Conventional Cutting	45	12	VIO
PGLA UND 4-0 X 75-19MM CU CT	PGN409570	4-0		19 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 4-0 X 75-26MM CU CT	PGN409580	4-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGLA UND 4-0 X 45-26MM CU CT	PGN409520	4-0		26 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA UND 4-0 X 45-19MM CU CT	PGN409510	4-0		19 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 4-0 X 45-19MM CU RC	PGN409386	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGLA UND 4-0 X 45-16MM CU RC	PGN405696	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	UND
PGLA VIO 4-0 X 75-16MM CU CT	PGN409443	4-0		16 mm 3/8 Circle Conventional Cutting	75	12	VIO

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)							
PGLA UND 4-0 X 75-17MM CU TC	PGN409825	4-0		17 mm 3/8 Circle Taper Cut	75	12	UND
PGLA VIO 4-0 X 90-45MM ST CT EX	PGN409711	4-0		45 mm Straight Conventional Cutting	90	12	VIO
PGLA VIO 4-0 X 70-16MM HC RB-3D	PGN402305Z	4-0		16 mm 1/2 Circle Round Body	70	36	VIO
PGLA UND 4-0 X 70-19MM CU RC-3D	PGN400422Z	4-0		19 mm 3/8 Circle Reverse Cutting	70	36	UND
PGLA UND 4-0 X 45-45MM ST CT	PGN409711USL	4-0		45 mm Straight Conventional Cutting	45	12	UND
PGLA VIO 4-0 X 45-13MM CU RC DN	PGN400559	4-0		13 mm 3/8 Circle Reverse Cutting Reverse Cutting	45	12	VIO
PGLA UND 4-0 X 45-13MM CU CT	PGN4083565U	4-0		13 mm 3/8 Circle Conventional Cutting	45	12	UND
PGLA VIO 4-0 X 75-30MM HC RB	PGN409430	4-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGLA UND 4-0 X 70-26MM HC RB	PGN400415U	4-0		26 mm 1/2 Circle Round Body	70	12	UND

Description	Code	Family
-------------	------	--------

MAXIMUS (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLACTIN 910)

PGLA VIO 4-0 X15 M CASSETTE SUTURE	PGNMC4015	Mitsu
PGLA VIO 3-0 X15 M CASSETTE SUTURE	PGNMC3015	Mitsu
PGLA VIO 2-0 X15 M CASSETTE SUTURE	PGNMC2015	Mitsu
PGLA VIO 0 X15 M CASSETTE SUTURE	PGNMC1015	Mitsu
PGLA VIO 1 X15 M CASSETTE SUTURE	PGNMC0115	Mitsu
PGLA VIO 2 X 25 M CASSETTE SUTURE	PGNMC0225	Mitsu
PGLA VIO 3 X25 M CASSETTE SUTURE	PGNMC0325	Mitsu
PGLA VIO 5 X25 M CASSETTE SUTURE	PGNMC0525	Mitsu
PGLA VIO 1 X 75 M CASSETTE SUTURE	PGNMC0175	Mitsu
PGLA VIO 1 X50 M CASSETTE SUTURE	PGNMC0150	Mitsu

MaxAbsorb CAS

Polyglycolic Acid Suture

Properties

- Retention of 70% - 80% knot tensile strength after 14 days
- Retention of 50% knot tensile strength after 21 days
- Complete mass absorption within 60 to 90 days

Indications

- Gastrointestinal procedures
- Gynaecology procedures
- Ophthalmic procedures
- Skin closure (intracutaneous, subcutaneous)
- Ligatures
- Abdominal Wall
- Muscles

Advantages

- High tensile strength
- Excellent knotting ability
- Smooth passage through tissue
- Easy handling





















Product Overview

Presentation/Format	: Needed suture packed in a Dozen box Non-needed suture packed in a cassette
Structure	: Multifilament, Braided
Colour	: Dyed, Undyed
Chemical composition	: Polyglycolic acid
Coating	: Polycaprolactone and Calcium Stearate
Sizes	: Violet - USP 8/0 (metric 0.4) - USP 2 (metric 5) : Undyed - USP 6/0 (metric 0.7) - USP 2 (metric 5)
Type of absorption	: Hydrolysis
Sterilisation	: Ethylene oxide

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXABSORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)

PGA VIO 2 X 75-30MM HC RB	MS029230	2		30 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2 X 90-40MM HC RB	MS2240N	2		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2 X 90-37MM HC RB	MS020348	2		37 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2 X 75-26MM HC RB	MS020238	2		26 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2 X 75-40MM HC RC	MS020481	2		40 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGA VIO 2 X 90-48MM HC RB	MS029252	2		48 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2 X 75-48MM HC RC	MS029297	2		48 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGA VIO 2 X 90-40MM HC RC	MS2478	2		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 2 X 75-37MM HC RB	MS020349	2		37 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2-0 X 90-30MM HC RB	MS2317	2-0		30 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2-0 X 140-30MM HC RB TC DN	MS2317DN	2-0		30 mm 1/2 Circle Round Body Taper Cut	140	12	VIO
PGA VIO 2-0 X 90-40MM HC RB	MS2345	2-0		40 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2-0 X 90-40MM HC RC	MS2382	2-0		40 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 2-0 X 75-40MM HC RB	MS2082544	2-0		40 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2-0 X 70-24MM CU RC	MS200443	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 2-0 X 90-30MM HC RC	MS200586	2-0		30 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 2 x 90 - 65MM HC RB	MS020880	2		65 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2 x 70 - 40MM CU RC	MS209241	2-0		40 mm 3/8 Circle Reverse Cutting	70	12	VIO
























Description	Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	-------------------------------------	--------------------	-----------	---------------

MAXABSORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)

PGA VIO 2 x 70 - 30MM CU RC	MS209242	2-0		30 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 2 x 70 - 48MM HC RC	MS209243	2-0		48 mm 1/2 Circle Reverse Cutting	70	12	VIO
PGA VIO 2 x 75 - 30MM HC RB EX	MS2317SL	2-0		30 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2 x 70 - 48MM HC RB	MS208060	2-0		48 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 2-0 X 100-26MM HC RB	MS209121	2-0		26 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 2-0 X 75-35MM HC RB	MS200737	2-0		35 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2-0 X 75-36MM HC TC	MS209363	2-0		36 mm 1/2 Circle Taper Cut	75	12	VIO
PGA VIO 2-0 X 90-36MM HC RB	MS209440	2-0		36 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 2-0 X 70-21MM HC RB	MS202824	2-0		21 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 2-0 X 75-36MM HC RB	MS209140	2-0		36 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 2-0 X 70-26MM HC RB	MS209123	2-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 2-0 X 100-36MM HC RB	MS209140XL	2-0		36 mm 1/2 Circle Round Body	100	12	VIO
PGA VIO 2-0 X 90-26MM HC RC	MS2869	2-0		26 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 2-0 X 70-17MM HC RB	MS200306	2-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 2-0 X 90-36MM HC RC	MS200466	2-0		36 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 2-0 X 70-24MM CU RB	MS202224	2-0		24 mm 3/8 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 70-20MM HC RB	MS2437	3-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 75-26MM HC CT	MS300774	3-0		26 mm 1/2 Circle Conventional Cutting	75	12	VIO
PGA VIO 3-0 X 70-26MM HC RB	MS300316	3-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 70-17MM HC RB	MS300667	3-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 75-26MM HC RB	MS300242	3-0		26 mm 1/2 Circle Round Body	75	12	VIO


















Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXABSORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)

PGA VIO 3-0 X 70-24MM CU RC	MS300442	3-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 3-0 X 70-30MM HC RB	MS309330	3-0		30 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 90-35MM HC RB	MS302738	3-0		35 mm 1/2 Circle Round Body	90	12	VIO
PGA VIO 3-0 X 70-24MM CU RB	MS302324	3-0		24 mm 3/8 Circle Round Body	70	12	VIO
PGA VIO 3-0 X 70-19MM CU RC	MS300398	3-0		19 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 3-0 X 70-22MM CU RC	MS302322	3-0		22 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 3-0 X 70-28MM CU RC	MS302328	3-0		28 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 3-0 X 75-24MM CU RC	MS2324N	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 3-0 X 90-26MM HC TC	MS302326	3-0		26 mm 1/2 Circle Taper Cut	90	12	VIO
PGA VIO 3-0 X 75-26MM CU RC	MS302326SL	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 3-0 X 75-30MM CU RC	MS302330	3-0		30 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 3-0 X 75-36MM HC RB	MS302736	3-0		36 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 3-0 X 90-22MM HC CT	MS2472	3-0		22 mm 1/2 Circle Conventional Cutting	90	12	VIO
PGA UND 3-0 X 70-19MM CU CT	MS309235	3-0		19 mm 3/8 Circle Conventional Cutting	70	12	UND
PGA UND 3-0 X 75-26MM CU RC	MS9932U	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	UND
PGA VIO 3-0 X 90-26MM HC RC	MS2868	3-0		26 mm 1/2 Circle Reverse Cutting	90	12	VIO
PGA VIO 3-0 X 75-22MM HC RB	MS309322	3-0		22 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 3-0 X 75-36MM HC RC	MS302336	3-0		36 mm 1/2 Circle Reverse Cutting	75	12	VIO
PGA VIO 4-0 X 70-20MM HC RB	MS2304	4-0		20 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 70-26MM HC RB	MS400315	4-0		26 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 45-19MM CU RC	MS400392	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA VIO 4-0 X 70-24MM CU RC	MS402672	4-0		24 mm 3/8 Circle Reverse Cutting	70	12	VIO
PGA VIO 4-0 X 75-19MM CU RC	MS400392XL	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	VIO

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXABSORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)

PGA VIO 4-0 X 70-24MM CU RB	MS402424	4-0		24 mm 3/8 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 70-17MM HC RB	MS400214	4-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 70 - 24MM CU RB	MS402424	4-0		24 mm 3/8 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 70-23MM CU RB	MS402423	4-0		23 mm 3/8 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 75-22MM HC RB	MS409422	4-0		22 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 4-0 X 70-16MM HC RB	MS2305	4-0		16 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 4-0 X 45-12MM CU RC	MS2494	4-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA UND 4-0 X 75-26MM CU CT	MS9918U	4-0		26 mm 3/8 Circle Conventional Cutting	75	12	UND
PGA VIO 5-0 X 45-16MM HC RB	MS2303	5-0		16 mm 1/2 Circle Round Body	45	12	VIO
PGA VIO 5-0 X 70-17MM HC RB	MS500303	5-0		17 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 70-13MM HC RB	MS502433	5-0		13 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 70-36MM HC RB	MS502536	5-0		36 mm 1/2 Circle Round Body	70	12	VIO
PGA VIO 5-0 X 45-13MM CU RB	MS502513	5-0		13 mm 3/8 Circle Round Body	45	12	VIO
PGA VIO 5-0 X 75-16MM HC RB	MS2303XL	5-0		16 mm 1/2 Circle Round Body	75	12	VIO
PGA VIO 5-0 X 75-16MM CU RC	MS500500XL	5-0		16 mm 3/8 Circle Reverse Cutting	75	12	VIO
PGA VIO 5-0 X 45-12MM CU RC	MS2493	5-0		12 mm 3/8 Circle Reverse Cutting	45	12	VIO
PGA VIO 5-0 X 45-15MM HC CT	MS500557	5-0		15 mm 1/2 Circle Conventional Cutting	45	12	VIO

Description	Code	Family
-------------	------	--------

MAXABSORB (SYNTHETIC ABSORBABLE BRAIDED COATED POLYGLYCOLIC ACID)

PGA VIO 2-0 X 15 M CASSETTE SUTURE	PGAMC2015	Megasorb
PGA VIO 0 X 15 M CASSETTE SUTURE	PGAMC1015	Megasorb
PGA VIO 1 X 15 M CASSETTE SUTURE	PGAMC0115	Megasorb

MaxiMus Plain CAS

Plain Catgut Suture

Properties

The tensile strength is maintained for not more than 10 to 12 days post implantation, and its absorption is complete within 90 to 100 days.

Indications

MERISOFT PLAIN™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability
- Rigidly controlled spinning & drying to ensure higher tensile strength



Product Overview

Primary packaging	: Packed in an aluminum foil pouch Non-needled suture packed in a cassette
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Yellow-tan
Sterilization	: Gamma radiation

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MAXIMUS PLAIN (PLAIN CATGUT SUTURE)							
PLAIN CATGUT 2-0 X 76-30MM CU RB	ME5033	2-0		30 mm 3/8 Circle Round Body	76	12	PLAIN
PLAIN CATGUT 2-0 X 100-37MM HC RB	ME204237	2-0		37 mm 1/2 Circle Round Body	100	12	PLAIN
PLAIN CATGUT 2-0 X 90-22MM HC RB	ME204222	2-0		22 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 2-0 X 90-26MM HC RB	ME204227	2-0		26 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 2-0 X 70-24MM CU RC	ME204224	2-0		24 mm 3/8 Circle Reverse Cutting	70	12	PLAIN
PLAIN CATGUT 2-0 X 90-50MM HC RB	ME207333	2-0		50 mm 1/2 Circle Round Body	50	12	PLAIN
PLAIN CATGUT 3-0 X 90-22MM HC RB	ME304322	3-0		22 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 3-0 X 90-17MM HC RB	ME304317	3-0		17 mm 1/2 Circle Round Body	90	12	PLAIN
PLAIN CATGUT 3-0 X 76-16MM CU CT	ME4285	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	PLAIN
PLAIN CATGUT 3-0 X 75-26MM HC RB	ME304327SL	3-0		26 mm 1/2 Circle Round Body	75	12	PLAIN
PLAIN CATGUT 3-0 X 76-26MM HC RB	ME4236	3-0		26 mm 1/2 Circle Round Body	76	12	PLAIN
PLAIN CATGUT 4-0 X 76-16MM CU RC	ME4283	4-0		16 mm 3/8 Circle Reverse Cutting	76	12	PLAIN
PLAIN CATGUT 2-0 X 152	SS2003	2-0		Non needed	152	12	PLAIN
PLAIN CATGUT 3-0 X 152	SS2002	3-0		Non needed	152	12	PLAIN

Description

Code

Family

MAXIMUS PLAIN (PLAIN CATGUT SUTURE)

PLAIN CATGUT 3-0 X100 M CASSETTE SUTURE	PLGMC30100	Merisoft Plain
PLAIN CATGUT 2-0 X100 M CASSETTE SUTURE	PLGMC20100	Merisoft Plain
PLAIN CATGUT 0 X100 M CASSETTE SUTURE	PLGMC10100	Merisoft Plain
PLAIN CATGUT 1 X75 M CASSETTE SUTURE	PLGMC0175	Merisoft Plain
PLAIN CATGUT 2 X50 M CASSETTE SUTURE	PLGMC0250	Merisoft Plain
PLAIN CATGUT 3 X50 M CASSETTE SUTURE	PLGMC0350	Merisoft Plain
PLAIN CATGUT 4 X30 M CASSETTE SUTURE	PLGMC0430	Merisoft Plain
PLAIN CATGUT 5 X30 M CASSETTE SUTURE	PLGMC0530	Merisoft Plain
PLAIN CATGUT 6 X30 M CASSETTE SUTURE	PLGMC0630	Merisoft Plain

MaxiMus Chromic CAS

Chromic Catgut Suture

Properties

MERISOFT CHROMIC™ gut suture is treated with a chromium salt solution to resist the enzymatic action, tensile strength retains until 10 - 14 days and absorption between 90-100 days.

Indications


MERISOFT CHROMIC™ suture is indicated for use in general soft tissue approximation and/or ligation.

Advantages

- Highest level of collagen guarantees predictable absorption
- Accurate polishing ensures smooth tissue passage
- Packed in liquid to reduce memory & increase pliability




















Product Overview

Primary packaging	: Packed in an aluminum foil pouch Non-needled suture packed in a cassette
Secondary packaging	: Primary pack is loaded in a secondary pouch made of medical grade paper & laminated polyethylene film
Sizes	: USP 6/0 (metric 1) - USP 2 (metric 6)
Colour	: Brown dyed 
Sterilization	: Gamma radiation

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMUS CHROMIC (CHROMIC CATGUT SUTURE)

CHROMIC CATGUT 2 X 100-45MM HC RB	ME4228	2		45 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-50MM HC RB	ME024250	2		50 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-40MM HC TRO	ME4244	2		40 mm 1/2 Circle Trocar	100	12	CHROMIC
CHROMIC CATGUT 2 X 100-65MM HC RB	ME4278	2		65 mm 1/2 Circle Round Body	100	12	CHROMIC
CHROMIC CATGUT 2 X 75-40MM HC RC	ME024244	2		40 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM CU RB	ME4216	2-0		30 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM HC RB	ME4241	2-0		30 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 152-36MM HC RCRB DN	ME4277	2-0		36 mm 1/2 Circle Reverse Cutting Round Body	152	12	CHROMIC
CHROMIC CATGUT 2-0 X 75-40MM HC RB	ME200991	2-0		40 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-26MM HC RB	ME204225	2-0		26 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-40MM HC RB	ME4245	2-0		40 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-45MM HC RB	ME4249	2-0		45 mm 1/2 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-30MM HC CT	ME4255	2-0		30 mm 1/2 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 76-50MM ST RB	ME4208	2-0		50 mm Straight Round Body	76	12	CHROMIC
CHROMIC CATGUT 2-0 X 75-45MM HC CT	ME200863	2-0		45 mm 1/2 Circle Conventional Cutting	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-22MM CU CT	ME4201	3-0		22 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-20MM HC RB	ME4237	3-0		20 mm 1/2 Circle Round Body	76	12	CHROMIC

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
CHROMIC CATGUT 3-0 X 75-26MM HC RB	ME300882	3-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-30MM HC RB	ME300172	3-0		30 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM CU RC	ME304326	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 3-0 X 38-16MM FE CT	ME4265	3-0		16 mm 5/8 Circle Conventional Cutting	38	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-16MM CU CT	ME4274	3-0		16 mm 3/8 Circle Conventional Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-12MM HC RC	ME4286	3-0		12 mm 1/2 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 76-26MM CU RC	ME4828	3-0		26 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 3-0 X 75-26MM HC RC	ME300886	3-0		26 mm 1/2 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-16MM CU RB	ME4048	4-0		16 mm 3/8 Circle Round Body	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-16MM CU RC	ME4287	4-0		16 mm 3/8 Circle Reverse Cutting	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 76-17MM CU TC	ME5047	4-0		17 mm 3/8 Circle Taper Cut	76	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-20MM HC RB	ME400181	4-0		20 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-17MM HC RB	ME400117	4-0		17 mm 1/2 Circle Round Body	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-19MM CU RC	ME400419	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	CHROMIC
CHROMIC CATGUT 4-0 X 38-12MM CU RC	ME4266	4-0		12 mm 3/8 Circle Reverse Cutting	38	12	CHROMIC
CHROMIC CATGUT 4-0 X 75-26MM HC RB	ME400121	4-0		26 mm 1/2 Circle Round Body	75	12	CHROMIC

Description

Code

Family

MAXIMUS CHROMIC (CHROMIC CATGUT SUTURE)

CHROMIC CATGUT 3-0 X100 M CASSETTE SUTURE	CHGMC30100	MeriSoft Chromic
CHROMIC CATGUT 2-0 X100M CASSETTE SUTURE	CHGMC20100	MeriSoft Chromic
CHROMIC CATGUT 0 X100M CASSETTE SUTURE	CHGMC10100	MeriSoft Chromic
CHROMIC CATGUT 1 X75 M CASSETTE SUTURE	CHGMC0175	MeriSoft Chromic
CHROMIC CATGUT 2 X 50 M CASSETTE SUTURE	CHGMC0250	MeriSoft Chromic
CHROMIC CATGUT 3 X50 M CASSETTE SUTURE	CHGMC0350	MeriSoft Chromic
CHROMIC CATGUT 4 X 30 M CASSETTE SUTURE	CHGMC0430	MeriSoft Chromic
CHROMIC CATGUT 5-0 X100M CASSETTE SUTURE	CHGMC50100	MeriSoft Chromic
CHROMIC CATGUT 4-0 X100M CASSETTE SUTURE	CHGMC40100	MeriSoft Chromic
CHROMIC CATGUT 4-0 X50 M CASSETTE SUTURE	CHGMC4050	MeriSoft Chromic
CHROMIC CATGUT 3-0 X50 M CASSETTE SUTURE	CHGMC3050	MeriSoft Chromic
CHROMIC CATGUT 2-0 X50 M CASSETTE SUTURE	CHGMC2050	MeriSoft Chromic
CHROMIC CATGUT 0 X75 M CASSETTE SUTURE	CHGMC1075	MeriSoft Chromic
CHROMIC CATGUT 0 X50 M CASSETTE SUTURE	CHGMC1050	MeriSoft Chromic
CHROMIC CATGUT 1 X50 M CASSETTE SUTURE	CHGMC0150	MeriSoft Chromic
CHROMIC CATGUT 4 X 25 M CASSETTE SUTURE	CHGMC0425	MeriSoft Chromic

MaxiMide CAS

Polyamide Suture

Indications


- Skin closure (intracutaneous, subcutaneous)
- Vascular procedures
- Ophthalmic procedures
- Neurosurgery

Advantages

- Flexible, easy to handle and tie
- Smooth passage through tissue
- Excellent histocompatibility




Product Overview

Presentation/Format	: Needed suture packed in a Dozen box Non-needed suture packed in a cassette
Structure	: Monofilament
Colour	: Black 
Chemical composition	: Polyamide 6 - 6.6
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Gamma radiation

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMIDE (POLYAMIDE SUTURE)

POLYAMIDE BLK 2 X 90-48MM HC RC	NYL023259	2		48 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2 X 75-45MM HC RB	NYL023245	2		45 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 2 X 90-90MM CU RC	NYL023290	2		90 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2 X 75-40MM HC RC	NYL022407	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 2 X 70-90MM CU CT	NYL02 3398	2		90 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 2 X 75-35MM HC RB	NYL023235	2		35 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 2 X 45-78MM CU CT	NYL020290	2		78 mm 3/8 Circle Conventional Cutting	45	12	BLACK
POLYAMIDE BU 2 X 100-90MM CU CT	NYL020798BU	2		90 mm 3/8 Circle Conventional Cutting	100	12	BLUE
POLYAMIDE BLK 2 X 75-40MM HC RC EX	NYL02 3240	2		40 mm 1/2 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 2-0 X 90-38MM CU RC	NYL202422	2-0		38 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-25MM CU RC	NYL206641	2-0		25 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-36MM HC RC	NYL203256	2-0		36 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 70-45MM CU RC	NYL20 3336	2-0		45 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 2-0 X 35-45MM CU RC	NYL20 3336SL	2-0		45 mm 3/8 Circle Reverse Cutting	35	12	BLACK
POLYAMIDE BLK 2-0 X 100-31MM HC RB	NYL200568	2-0		31 mm 1/2 Circle Round Body	100	12	BLACK
POLYAMIDE BLK 2-0 X 45-26MM CU RC	NYL200664	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 2-0 X 90-26MM CU RC	NYL201626	2-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-30MM CU RC	NYL203232	2-0		30 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 75-24MM CU CT	NYL203224	2-0		24 mm 3/8 Circle Conventional Cutting	75	12	BLACK

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------








MAXIMIDE (POLYAMIDE SUTURE)

POLYAMIDE BLK 2-0 X 90-26MM HC RC	NYL209026	2-0		26 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 70-30MM CU RB	NYL207030	2-0		30 mm 3/8 Circle Round Body	70	12	BLACK
POLYAMIDE BLK 2-0 X 90-55MM ST CT	NYL200776	2-0		55 mm Straight Conventional Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-30MM HC RC	NYL203682	2-0		30 mm 1/2 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 2-0 X 90-60MM ST CT	NYL200628	2-0		60 mm Straight Conventional Cutting	90	12	Black
POLYAMIDE BLK 2-0 X 75-19MM CU RC	NYL2033210	2-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 2-0 X 45-24MM CU RC	NYL2093368	2-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 90-25MM CU RB	NYL306631	3-0		25 mm 3/8 Circle Round Body	90	12	BLACK
POLYAMIDE BLK 3-0 X 90-35MM CU RC	NYL301627	3-0		35 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 70-26MM CU RC	NYL30 3328	3-0		26 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 35-26MM CU RC	NYL30 3328SL	3-0		26 mm 3/8 Circle Reverse Cutting	35	12	BLACK
POLYAMIDE BLK 3-0 X 70-16MM CU CT	NYL30 3321	3-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 3-0 X 45-19MM CU RC	NYL301621	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 90-26MM CU RC	NYL301625	3-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 90-19MM CU RC	NYL301622	3-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 3-0 X 75-35MM HC RB	NYL303335	3-0		35 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-25MM CU CT	NYL303326	3-0		25 mm 3/8 Circle Conventional Cutting	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-40MM HC RB	NYL305340	3-0		40 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 3-0 X 70-60MM ST CT	NYL30 3388	3-0		60 mm Straight Conventional Cutting	70	12	BLACK
POLYAMIDE BU 3-0 X 45-19MM CU RC	NYL3013498	3-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 3-0 X 75-26MM CU RC	NYL301685	3-0		26 mm 3/8 Circle Reverse Cutting	75	12	BLUE

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXIMIDE (POLYAMIDE SUTURE)

POLYAMIDE BU 3-0 X 75-60MM ST CT	NYL306441	3-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYAMIDE BU 3-0 X 45-26MM CU RC	NYL301625BU	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 3-0 X 45-24MM CU RC	NYL301857	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BLK 3-0 X 100-30MM CU RC-3D	NYL303330Z	3-0		30 mm 3/8 Circle Reverse Cutting	100	36	BLACK
POLYAMIDE BLK 3-0 X 90-22MM HC RB	NYL309022	3-0		22 mm 1/2 Circle Round Body	90	12	BLACK
POLYAMIDE BLK 3-0 X 45-26MM CU RC	NYL301625SL	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 3-0 X 75-24MM CU RC	NYL305669	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 3-0 X 75-60MM ST CT	NYL300627	3-0		60 mm Straight Conventional Cutting	75	12	Black
POLYAMIDE BU 3-0 X 100-26MM CU CT	NYL3013500	3-0		26 mm 3/8 Circle Conventional Cutting	100	12	Blue
POLYAMIDE BLK 3-0 X 100-30MM CU RC	NYL303330	3-0		30 mm 3/8 Circle Reverse Cutting	100	12	Black
POLYAMIDE BLK 3-0 X 75-20MM HC RC	NYL303370	3-0		20 mm 1/2 Circle Reverse Cutting	75	12	Black
POLYAMIDE BLK 4-0 X 90-16MM CU CT	NYL40 3319	4-0		16 mm 3/8 Circle Conventional Cutting	90	12	BLACK
POLYAMIDE BLK 4-0 X 45-11MM CU RC	NYL401612	4-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU RC	NYL401619	4-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLACK
POLYAMIDE BLK 4-0 X 45-19MM CU RC	NYL401620	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 45-16MM CU RC	NYL401861	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLACK
POLYAMIDE BLK 4-0 X 90-19MM CU RC	NYL401618	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 4-0 X 75-19MM CU CT	NYL403419	4-0		19 mm 3/8 Circle Conventional Cutting	75	12	BLACK
POLYAMIDE BLK 4-0 X 75-28MM HC RB	NYL405428	4-0		28 mm 1/2 Circle Round Body	75	12	BLACK
POLYAMIDE BLK 4-0 X 60-19MM CU CT	NYL40 3357	4-0		19 mm 3/8 Circle Conventional Cutting	60	12	BLACK
POLYAMIDE BLK 4-0 X 38-10MM CU RC	NYL40 3326	4-0		10 mm 3/8 Circle Reverse Cutting	38	12	BLACK
POLYAMIDE BLK 4-0 X 70-19MM CU CT	NYL40 3567	4-0		19 mm 3/8 Circle Conventional Cutting	70	12	BLACK

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MAXIMIDE (POLYAMIDE SUTURE)							
POLYAMIDE BU 4-0 X 45-11MM CU RC	NYL401612BU	4-0		11 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BU 4-0 X 45-19MM CU RC	NYL401620BU	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYAMIDE BLK 4-0 X 70-16MM CU CT	NYL40 3318	4-0		16 mm 3/8 Circle Conventional Cutting	70	12	BLACK
POLYAMIDE BLK 4-0 X 70-10MM CU RC	NYL40 3326XL	4-0		10 mm 3/8 Circle Reverse Cutting	70	12	BLACK
POLYAMIDE BLK 4-0 X 45-19MM CU RC-3D	NYL401620Z	4-0		19 mm 3/8 Circle Reverse Cutting	45	36	BLACK
POLYAMIDE BLK 4-0 X 90-24MM CU RC	NYL403424XL	4-0		24 mm 3/8 Circle Reverse Cutting	90	12	BLACK
POLYAMIDE BLK 4-0 X 75-18MM CU RB	NYL403418	4-0		18 mm 3/8 Circle Round Body	75	12	BLACK

Description	Code	Family
MAXIMIDE (POLYAMIDE SUTURE)		
NYLON BLACK 3-0 X100 MCASSETTE SUTURE	NYLMC30100	Filamide
NYLON BLACK 2-0 X100 MCASSETTE SUTURE	NYLMC20100	Filamide

MaxProp CAS

Polypropylene Suture

Indications

- Skin closure (intracutaneous, subcutaneous)
- Cardiovascular procedures
- Neurological procedures
- Ophthalmic procedures

Advantages

- Smooth passage through tissue
- Minimal acute inflammatory reaction in tissues
- Unsurpassed strength
- Optimal elasticity and elongation properties
- Available in Less Memory Tray (LMT)





















Product Overview

Presentation/Format	: Needed suture packed in a Dozen box Non-needed suture packed in a cassette
Structure	: Monofilament
Colour	: Blue
Chemical composition	: Polypropylene
Sizes	: USP 10/0 (metric 0.2) - USP 2 (metric 5)
Sterilization	: Ethylene oxide

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------



















MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 2 X 75-45MM HC RC	PPL028245	2		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2 X 75-45MM HC RB	PPL028275	2		45 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM CU RC	PPL208689	2-0		40 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-90MM CU RC	PPL200987	2-0		90 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC TC	PPL208852	2-0		40 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB	PPL20 844	2-0		30 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC TC DN	PPL20 8977	2-0		26 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-17.5MM HCTC DN	PPL20 8937	2-0		17.5 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC RB	PPL208523	2-0		26 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 75-26MM HC RB	PPL208833	2-0		26 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-24MM CU RC	PPL208685	2-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-26MM HC CT	PPL208978	2-0		26 mm 1/2 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-60MM ST CT	PPL200623	2-0		60 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-24MM CU RC	PPL208685XL	2-0		24 mm 3/8 Circle Reverse Cutting	100	12	BLUE
POLYPROPYLENE BU 2-0 X 75-26MM HC RB DN	PPL208834	2-0		26 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-36MM HC RB	PPL208236	2-0		36 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-70MM ST CT	PPL200753	2-0		70 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 75-45MM HC RC	PPL208245	2-0		45 mm 1/2 Circle Reverse Cutting	75	12	BLUE

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------



















MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 2-0 X 75-45MM HC RB	PPL208275	2-0		45 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 70-30MM CU RB	PPL208230	2-0		30 mm 3/8 Circle Round Body	70	12	BLUE
POLYPROPYLENE BU 2-0 X 70-25MM HC TC	PPL20 807	2-0		25 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 2-0 X 75-20MM HC RB	PPL2085200	2-0		20 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-30MM CU RC	PPL207697	2-0		30 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-60MM ST CT	PPL200631	2-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-31MM HC RB	PPL208231	2-0		31 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 75-30MM HC RB	PPL20844SL	2-0		30 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-31MM HC RB DN	PPL208526	2-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-31MM HC TC DN	PPL208850	2-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 45-26MM CU RC	PPL208026	2-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC TC DN	PPL208852DN	2-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-75MM ST CT	PPL200752	2-0		75 mm Straight Conventional Cutting	100	12	BLUE
POLYPROPYLENE BU 2-0 X 90-55MM ST CT	PPL208571	2-0		55 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-40MM HC RB DN	PPL20 2040DN	2-0		40 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 2-0 X 70-30MM HC RB	PPL20 841	2-0		30 mm 1/2 Circle Round Body	70	12	BLUE

Description	Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
MAXPROP (POLYPROPYLENE SUTURE)						
POLYPROPYLENE BU 2-0 X 90-90MM CU CT	PPL200988	2-0		90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-26MM HC TC	PPL208926	2-0		90	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RC	PPL208931	2-0		90	12	BLUE
POLYPROPYLENE BU 2-0 X 75-70MM ST RB DN	PPL208400	2-0		75	12	BLUE
POLYPROPYLENE BU 2-0 X 90-30MM HC RB DN	PPL20844DN	2-0		90	12	BLUE
POLYPROPYLENE BU 2-0 X 100-40MM CU CT	PPL208690	2-0		100	12	BLUE
POLYPROPYLENE BU 3-0 X 70-22MM CU CT	PPL30 018	3-0		70	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM CU CT	PPL30 800	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM HC RB	PPL30 825	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 70-25MM HC RB	PPL30 838	3-0		70	12	BLUE
POLYPROPYLENE BU 3-0 X 70-20MM HC RB	PPL30 887	3-0		70	12	BLUE
POLYPROPYLENE BU 3-0 X 70-60MM ST CT	PPL30 888	3-0		70	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC RB DN	PPL30 8522	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17MM HC RB	PPL308558	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM HC TC DN	PPL308667	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-22MM HC RB	PPL308770	3-0		75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC TC DN	PPL308976	3-0		90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-19MM CU RC	PPL308665	3-0		45	12	BLUE

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 3-0 X 45-24MM CU RC	PPL308684	3-0		24 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 75-26MM HC RB	PPL308832	3-0		26 mm 1/2 Circle Round Body	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-19MM CU RC	PPL308687	3-0		19 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM HC RB	PPL308559	3-0		16 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-25MM HC TC	PPL308977	3-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-30MM HC RB	PPL308525	3-0		30 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-16MM CU RC	PPL308549	3-0		16 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM CU RC	PPL308685	3-0		26 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-17MM HC TC	PPL308317	3-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB	PPL308336	3-0		36 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 75-22MM HC RB DN	PPL308771	3-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-60MM ST CT	PPL3013413	3-0		60 mm Straight Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 75-16MM CU CT	PPL3013453	3-0		16 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC RB	PPL308534	3-0		26 mm 1/2 Circle Round Body	120	12	BLUE
POLYPROPYLENE BU 3-0 X 75-24MM CU RC	PPL307694	3-0		24 mm 3/8 Circle Reverse Cutting	75	12	BLUE
POLYPROPYLENE BU 3-0 X 90-36MM HC RB DN	PPL307499	3-0		36 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 120-26MM HC TC DN	PPL307778	3-0		26 mm 1/2 Circle Taper Cut Taper Cut	120	12	BLUE
POLYPROPYLENE BU 3-0 X 45-26MM CU RC	PPL308025	3-0		26 mm 3/8 Circle Reverse Cutting	45	12	BLUE

Description	Code	Size	Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	-------------------------------------	--------------------	-----------	---------------

MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 3-0 X 90-31MM HC RB DN	PPL308525DN	3-0		31 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-22MM HC CT	PPL308731	3-0		22 mm 1/2 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 45-19MM CU CT	PPL308878	3-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HC RB DNV	PPL308355	3-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-31MM HC TC DN	PPL308849DN	3-0		31 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-40MM HC TC DN	PPL308851	3-0		40 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-26MM HCRB TRODN	PPL308833	3-0		26 mm 1/2 Circle Round Body TRO	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-20MM HC RBTC DN	PPL309557	3-0		20 mm 1/2 Circle Round Body Taper Cut	90	12	BLUE
POLYPROPYLENE BU 3-0 X 90-60MM ST CT	PPL300621	3-0		60 mm Straight Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 3-0 X 45-26MM CU CT EX	PPL308021	3-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 3-0 X 75-36MM CU CT	PPL300626	3-0		36 mm 3/8 Circle Conventional Cutting	75	12	BLUE
POLYPROPYLENE BU 4-0 X 90-19MM CU RC	PPL408015	4-0		19 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-25MM CU CT	PPL408879	4-0		25 mm 3/8 Circle Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB DN	PPL40 8521	4-0		26 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB	PPL408355	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 4-0 X 90-17MM HC RB	PPL408558	4-0		17 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC RB	PPL408761	4-0		20 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC CT DN	PPL408840	4-0		20 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC CT DN	PPL408846	4-0		26 mm 1/2 Circle Conventional Cutting Conventional Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 45-19MM CU RC	PPL408682	4-0		19 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-20MM CU RC	PPL408683	4-0		20 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-18MM CU RB DN	PPL408861	4-0		18 mm 3/8 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-24MM CU RC	PPL408684	4-0		24 mm 3/8 Circle Reverse Cutting	90	12	BLUE
POLYPROPYLENE BU 4-0 X 70-19MM CU RC	PPL408682XL	4-0		19 mm 3/8 Circle Reverse Cutting	70	12	BLUE
POLYPROPYLENE BU 4-0 X 90-25MM HC TC	PPL408425	4-0		25 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 75-22MM HC RB DN	PPL40 8526DN	4-0		22 mm 1/2 Circle Round Body Round Body	75	12	BLUE
POLYPROPYLENE BU 4-0 X 70-17MM HC TC	PPL40 835	4-0		17 mm 1/2 Circle Taper Cut	70	12	BLUE
POLYPROPYLENE BU 4-0 X 45-16MM CU RC	PPL408884	4-0		16 mm 3/8 Circle Reverse Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-26MM CU CT	PPL408020SL	4-0		26 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 45-19MM CU CT	PPL408011	4-0		19 mm 3/8 Circle Conventional Cutting	45	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC TC	PPL40 857	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-17MM HC TC EX	PPL400835	4-0		17 mm 1/2 Circle Taper Cut	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-26MM HC RB EX	PPL408831	4-0		26 mm 1/2 Circle Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC TC DN	PPL408841DN	4-0		20 mm 1/2 Circle Taper Cut Taper Cut	90	12	BLUE

Description	Code	Size		Needle Profile Shape, Type & Length	Suture Length (cm)	Pack Size	Suture Colour
-------------	------	------	--	-------------------------------------	--------------------	-----------	---------------

MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 4-0 X 45-20MM CU RC-3D	PPL408683Z	4-0		20 mm 3/8 Circle Reverse Cutting	45	36	BLUE
POLYPROPYLENE BU 4-0 X90-26MM HCRB EX-3D	PPL408831Z	4-0		26 mm 1/2 Circle Round Body	90	36	BLUE
POLYPROPYLENE BU 4-0 X 90-20MM HC RB DN	PPL408761DN	4-0		20 mm 1/2 Circle Round Body Round Body	90	12	BLUE
POLYPROPYLENE BU 4-0 X 90-22MM HC RB	PPL408762	4-0		22 mm 1/2 Circle Round Body	90	12	BLUE

Description	Code	Family
-------------	------	--------

MAXPROP (POLYPROPYLENE SUTURE)

POLYPROPYLENE BU 1 X 25 M CASSETTE SUTUR	PPLMC0125	FILAPROP
POLYPROPYLENE BU 2 X 25 M CASSETTE SUTUR	PPLMC0225	FILAPROP
POLYPROPYLENE BU 0 X 25 M CASSETTE SUTUR	PPLMC1025	FILAPROP
POLYPROPYLENE BU 2-0 X 25 M CASSETTE SUT	PPLMC2025	FILAPROP
POLYPROPYLENE BU 3-0 X 25 M CASSETTE SUT	PPLMC3025	FILAPROP
POLYPROPYLENE BU 4-0 X 25 M CASSETTE SUT	PPLMC4025	FILAPROP